



Anodyne® Therapy Consumer System

Model 120

Manufacturer:
Anodyne Therapy, LLC
Address: 301B Mears Blvd.
Oldsmar, FL 34677
Tel: 1-800-521-6664
Fax: 1-800-835-4581
Website: www.anodynetherapy.com



Anodyne® Therapy and MIRE™
are registered trademarks of Anodyne Therapy, LLC.
US Patent No. 5,358,503. European Patent No. 0741594

© 2018 Anodyne® Therapy, LLC
All Rights Reserved

Important Safety Information and Instructions

Read Entire Booklet Before Operating



The Anodyne® Therapy System is a noninvasive medical device that delivers Monochromatic Infrared Energy/ (MIRE™) through infrared light-emitting diodes. These diodes are mounted in flexible Therapy Pads, and emit infrared light at a wavelength of 890 nm - increasing local circulation and reducing pain, stiffness and muscle spasm where applied.

Anodyne® Therapy Systems received clearance from the US FDA in 1994 and are compliant to the ISO 13485 International Standard.

Review this manual in its entirety prior to using this system and contact Anodyne Therapy, LLC at 1-800-521-6664 with any questions, comments or concerns.

INDEX

Cleaning, 15

Contraindications, 4

Control Unit, 5, 7-8, 16-19

Customer Service, 2, 4-5, 6, 8, 11, 17

Device Disposal, 5

Pain Patches, 4, 18

Repairs, 17

Return Policy, 22

Single Patient Use, 6, 9, 14

Straps, 5, 7, 9-10, 12, 15, 19

Therapy Pad Placement, 4, 9-14, 19

Therapy Pad Pressure, 4, 10-13, 19

Topical Agents, 4, 18

Troubleshooting, 17

Treatment Frequency, 12-13

Treatment Time, 12-13

Warnings, Precautions, 4-9, 11-13, 15-16, 18-19

Warranty, 16, 22



LIMITED WARRANTY

Anodyne Therapy, LLC (“Manufacturer”) warrants the Anodyne® Model 120 (“The Product”) to the **immediate purchaser only** as follows:

Limited Warranty

Manufacturer warrants that The Product sold hereunder will be free from defects in material and workmanship for a period of either: two (2) years (New Systems) or one (1) year (Refurbished Systems) from the date of purchase with normal use. If the defects are of such type and nature as to be covered by this warranty, **Manufacturer shall, at its option, either repair or replace the damaged product at its expense, (minus customer shipment fees of products to manufacturer).** This warranty does not cover any Products that have been abused, misused, or tampered with in any way. This warranty does not cover Products damaged in a natural disaster or flood; nor does it cover theft. **This limited warranty is not transferable, and only applies to purchases direct from Anodyne Therapy LLC, or its authorized distributors.**

THIS WARRANTY IS IN LIEU OF OTHER WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR USE, WHICH ARE HEREBY SPECIFICALLY DISCLAIMED.

This limited warranty gives you specific legal rights. You may have other rights, which vary from state to state. To the extent allowed by applicable law, in no event shall manufacturer be liable for any incidental, consequential, special, indirect, punitive or exemplary damages or lost profits from any breach of warranty.

TABLE OF CONTENTS

General Warnings and Precautions.....	4-6
Purpose of Device (Indications for Use)	6
Device and Parts Description.....	7
Environmental Conditions Affecting Use	7
Instructions for Use	8-10
Treatment Guidelines	11-14
Cleaning	15
Maintenance	16
Storage.....	16
Troubleshooting	17
Return a Unit For Repair.....	17
Summary	18-19
Declarations	20-21
Warranty and Return Policy.....	22
Index	23





CONTRAINDICATIONS

The Anodyne® Therapy System, Model 120 should not be used directly over or near the womb during **pregnancy** or directly over **an active cancer**. The use of this device with these types of conditions are unknown.

GENERAL WARNINGS AND PRECAUTIONS

WARNING: Use carefully. May cause serious burns. Do not use over sensitive areas or in the presence of poor circulation. The unattended use of the Model 120 by children or incapacitated persons may be dangerous.

To minimize the risks of burns, please read and follow these important precautions:

Therapy Pad Placement	<ul style="list-style-type: none"> Never use the Therapy Pads for more than 45 minutes in one area. If skin is very dry, thin or sensitive, use for the minimum times listed in the <i>Instructions for Use</i>. Slightly lift your Therapy Pads every 10 minutes to check for burns. If there is evidence of a burn, discontinue treatment and call 1-800-521-6664. Never sleep during an active treatment session. We recommend setting a timer to avoid over treatment.
Therapy Pad Pressure	<ul style="list-style-type: none"> Do not stand on or put excessive pressure on either of the Therapy Pads during treatment as this will increase the risk of burns. Use caution when treating over bony areas such as ankles, knees, elbows, and tops of feet and hands. These areas have less tissue and are more susceptible to burns.
Therapy Pad Warmth	<ul style="list-style-type: none"> It is normal to experience warmth during treatment. However, if any of the Therapy Pads feel too hot for comfort, discontinue use and call 1-800-521-6664.
Skin Color Changes	<ul style="list-style-type: none"> It is normal for the treated area to be slightly pink after treatment. However, if area is red like a sunburn and stays that way for more than 8 hours - do not treat again prior to contacting 1-800-521-6664.
Topical Agents	<ul style="list-style-type: none"> Do not use Therapy Pads over pain patches like Duragesic®, or residual of; or topical heating agents like BenGay® or Icy Hot®, as this may increase the likelihood of a burn.
Other	<ul style="list-style-type: none"> Do not saturate, soak, or immerse components in water or liquid.

DECLARATION - ELECTROMAGNETIC IMMUNITY

The Anodyne® Therapy System, Model 120, is intended for use in the electromagnetic environment specified below. The user of this product should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+ - 6 kV contact + - 8 kV air	+ - 6 kV contact + - 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-5	+ - 2 kV for power supply lines + - 1 kV for i/o lines	+ - 2 kV for power supply lines + - 1 kV for i/o lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ - 1 kV differential mode + - 2 kV common mode	+ - 1 kV differential mode + - 2 kV common mode	Main power should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% Ut (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) <5% UT (>95% dip in UT) for 5 s	<5% Ut (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) <5% UT (>95% dip in UT) for 5 s	Main power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field. ISO 61000-4-8	3 A/m	0.3 A/m	If image distortion occurs, it may be necessary to position Model 120 further from the sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

NOTE: UT is the A.C. main voltage prior to application of the test level.



DECLARATION - ELECTROMAGNETIC EMISSIONS

The Anodyne® Therapy System, Model 120, is intended for use in the electromagnetic environment specified below. The user of this product should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions	Group 1	The Model 120 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The Model 120 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supply buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage Fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

DECLARATION - PHOTO BIOLOGICAL EMISSIONS

Each Therapy Pad contains 60 - 890nm light emitting diodes (LEDs) that emit infrared energy at a radiant intensity of 18.4 mW/sr. Based upon the supplier's LED specifications, the photo biological emissions from the Anodyne Therapy System Model 120 meet the eye safety guidelines of the IEC60825-1- Class 1, IEC62471- Exempt, and the EU Directive 2006/25/EC.



GENERAL WARNINGS AND PRECAUTIONS

CAUTION: THIS IS AN ELECTRICAL DEVICE AND CAN CAUSE ELECTRICAL SHOCK.

Electrical Precautions

- The device is intended to be connected with a medical grade power supply plugged into a standard receptacle.
- This equipment is intended for use with Anodyne® Therapy pad cables only.
- The Anodyne® Therapy System needs to be installed and put into service according to the information supplied with this instruction manual.
- No modification of this equipment is allowed.
- Portable and mobile RF communications equipment can affect the proper operation of medical electronic equipment. Call 1-800-521-6664 before using this device if you, or any other member of your family, is using other medical electronic equipment in your home.
- Use only the supplied attachment straps, paper tape, gravity or applicators approved in this manual to maintain the Therapy Pad position. Do not use pins or other metallic materials as a means of attachment.
- Like any other electrical device, do not use this system in water, near water or while wet. Use of this system in these situations could cause electrical shock.
- Do not use a portable electric generator as a power source as this may damage the electrical components of your device.
- The Anodyne® Therapy System, Model 120 is an electronic device in which a portion of the electronics may be hazardous. As such, the device should be disposed of in an environmentally responsible manner. The device should be disposed of in a hazardous waste management facility in your local area at the end of its useful life.

Markings and Symbols

	Attention, consult accompanying documents
	Off (power: disconnection from the mains)
	On (power: connection to the mains)
	Direct current
	Alternating current
	Dispose at hazardous waste facility



GENERAL WARNINGS AND PRECAUTIONS

GENERAL WARNINGS

- **The Model 120 emits heat and may cause burns!**
If any of the Therapy Pads feel too hot, immediately stop treatment by depressing the On/Off switch to the off position, remove the Therapy Pads, and call 1-800-521-6664.
- **The Model 120 has been designed for single patient use only.** Sharing this device with others may result in the spread of bacteria, viruses and other micro organisms from one person to another, and may spread disease.
- If there is an open wound present in/near the treatment area, cover the Therapy Pads or the treatment area (entire foot if using the soft shoe) with a clear plastic barrier to keep this system clean. Discard plastic barrier after use.
- **If you are diabetic, as your pain decreases, your activity level may increase, and this could cause hypoglycemia.** Therefore, be sure to monitor your blood sugars regularly. Use of the Model 120 alone does not cause hypoglycemia.
- Do not leave your Model 120 unattended around small children or pets. This could lead to serious shock, cuts and or infections.
- Do not use if device materials cause skin irritation.

PURPOSE OF THE DEVICE (INDICATIONS FOR USE)


- The purpose of the Anodyne® Model 120 is for the temporary relief of pain, stiffness and muscle spasm and for temporarily increasing blood circulation.

SUMMARY


Basic Operating Instructions:

Set Up: Place the Control Unit on the floor next to you. Plug the Therapy Pad cords into the Connectors on the front of the Control Unit. Plug Power Supply into an electrical outlet. Place the Therapy Pads over the area to be treated and attach with supplied attachment straps.


If using the Lower Extremity Treatment Application System, one Therapy Pad should be positioned horizontally across the toes, and the second Therapy pad should be positioned lengthwise across the ball of the foot. Place two Therapy Pad on either side of your calf where your pain or poor circulation begins and secure with the Leg Attachment Strap.

 **Warning to decrease the risk of burns,** do not over tighten Attachment Straps, stand or put pressure on Foot Therapy Pad. Do not treat on the top of the foot or hand and be careful treating over bony areas.

Turn on Control unit to begin treatment and depress the On/Off Switch on the Control Unit to the On position. If the Green Indicator Light turns on the device is on and operational.

 **Warning:** Immediately remove Therapy Pads if they feel too warm for comfort as this may increase the risk of burns.

Remove Pads at the end of treatment: Loosen Attachment Straps of the Calf and Foot Therapy Pads and slip off of the leg and foot. Turn off Power Adapter, unplug from electrical outlet, and store in a safe place away from children and pets.

 **Warning:** Do not store in a wet or excessively damp location as this may cause rust or corrosion of the electrical components or wiring and may result in electrical shock or impair function of this device.



SUMMARY

This summary is not intended to replace the complete user instructions. The entire manual should be read and understood before operating the Model 120.



Warning and Cautions

This device emits heat and may cause burns. Risk of burns is reduced in the following ways:

- Never use pain patches or topical heating agents before, during, or after treatment.
- Stop treatment if the Control Unit or Therapy Pads feel too warm.
- Apply the Therapy Pads exactly as described throughout this manual.

This is an electrical device. To avoid electrical shock:

- Do not use in water, near water or while wet.
- Do not use if electrical cords are frayed or you can see exposed wires.
- Do not saturate, soak or immerse components in water or liquid.
- Do not store in a wet or excessively damp location-this may cause rust or corrosion of electrical components or wiring and result in electrical shock or impair proper function of device.

COMPONENTS

The Anodyne® Therapy System, Model 120 consists of the Main Control Unit, Therapy Pads, Attachment Straps, Velcro Strips, and an optional Lower Extremity Treatment System. Also included is a carrying case, which is not pictured.



NORMAL ENVIRONMENTAL CONDITIONS

- An ambient temperature range of +10 degrees C to +40 degrees C (50° F to 104° F).
- A relative humidity range of 30% to 75%.
- An atmospheric pressure range of 500 hPa to 1060 hPa.

ENVIRONMENTAL CONDITIONS FOR TRANSPORTATION AND STORAGE

- An ambient temperature range of -40 degrees C to 70 degrees C (-40° F to 158° F).
- A relative humidity range of 10 to 100RH including condensation. Be sure not to store in excessively damp location.
- An atmospheric pressure range of 500 hPa to 1060 hPa.
- Device is packaged to maintain cleanliness and minimize egress of water.



Do not operate in a wet environment as this may cause electrical shock.



INSTRUCTIONS FOR USE

SETTING UP YOUR SYSTEM



Warning: If electrical wires are exposed, discontinue system use and call Anodyne Therapy, LLC at 1-800-521-6664 to arrange for service. Exposed wires may be unsafe and cause electrical shock.

Step 1: Remove the device from the packaging and protective plastic bag. The plastic bag may be discarded or retained for product storage between uses.

Step 2: Read the entire user manual before using this device.

Step 3: If you are treating your feet and legs, remove socks or other hosiery - this will ensure you obtain the best results.

Step 4: Place the Control Unit on the floor or on a table next to either an easy chair or couch where you will treat yourself. Do not position the device where it could fall and cause injury. This device is portable and does not have to be set up in the same location each time.

Step 5: Plug the Power Supply into an active electrical outlet.



Step 6: Turn on the Control Unit by depressing the On/Off Switch to the On position. The Green Indicator Light will turn on (if not, please see Troubleshooting).

Step 7: The system is now ready for use. Information on pad placement and treatment times is found in the Instructions for Use on pages 12 and 13 of this user manual.

Step 8: Ending Your Treatment Session: After your treatment is completed, turn the On/Off switch to the OFF position, unplug the Power Supply from the wall, and store in a safe, dry and cool place away from children and pets until the next use.

TROUBLESHOOTING GUIDE

Problem	Troubleshooting
Only one green light illuminates on front of the Control Unit	<ul style="list-style-type: none"> Manipulate the Therapy Pad wires, does this cause the green lights to re-appear or flicker? Discontinue system use and call Anodyne Therapy, LLC at 1-800-521-6664 to arrange for service.
Green light on Power Supply does not illuminate/no power	<ul style="list-style-type: none"> Make sure that the Power Supply plug is securely inserted into an electric outlet. Remove any extension cords in use, and plug the system directly into the electrical outlet of the wall. Discontinue system use and call Anodyne Therapy, LLC at 1-800-521-6664 to arrange for service.
Therapy Pads are damaged	<ul style="list-style-type: none"> Wires exposed or Therapy Pad materials are separating. Discontinue system use and call Anodyne Therapy, LLC at 1-800-521-6664 to arrange for service.
Control Unit or Therapy Pads get too hot	<ul style="list-style-type: none"> A therapeutic warmth is normal, however, this should not cause discomfort. Have a friend or family member confirm heat. If problem persists, discontinue system use and call Anodyne Therapy, LLC at 1-800-521-6664 to arrange for service.

RETURNING A UNIT FOR REPAIR

To send a unit to Anodyne Therapy, LLC for repair, you must first call 1-800-521-6664 and obtain a return authorization (RA). The following information must be supplied to Anodyne Therapy, LLC before an RA number can be issued:

- 1) User name and address
- 2) Device serial number
- 3) A description of the problem with the unit.

Place return items in a closed plastic bag; place bagged items in a shipping box. DO NOT RETURN straps, velcro attachment tabs, soft shoes or carrying case.



MAINTENANCE

No special maintenance procedures are required, however, proper care of the Model 120 using the following guidelines will extend the useful life of this product.

1. **Be careful not to drop your Control Unit on a hard surface**, step on it or drop a heavy object on it - the as this could damage the internal circuitry.
2. **Do not remove** the system Power Supply from an electrical outlet by pulling the cords - but rather by pulling by the electrical plug.
3. **Do not remove** the system Power Supply from the back of the Control Unit, as these components are soldered together.
4. **Do not** wrap cords around the Therapy Pads, Control unit or Power Supply. This may cause the Therapy Pads to separate or cause the cords to fray.
5. **Do not** carry or pull the system by Therapy Pads or electrical cords. Doing so may cause components to separate and or fray.
6. When your treatment session is complete, **always turn your system off** using the On/Off Switch on the back of the Control Unit and unplug the Power Supply from the electric outlet.

STORAGE

When not in use, store your Model 120 in a cool, dry area away from access by either children or pets. Storage in a wet or excessively damp location may cause corrosion of electrical components or wiring that may result in electrical shock or impair the proper function of the device.

Otherwise, no special requirements exist for the transport and storage of the device. The device is shipped in a poly bag as a precaution from water damage, but this is not a requirement to transport the device.

INSTRUCTIONS FOR USE

SETTING UP YOUR SOFT SHOE

Never stand or walk in the Soft Shoe as this may increase your risk for a burn.



The Soft Shoe is not meant to be used with multiple patients. Do not share your Soft Shoe with others.


This system comes with the Therapy Pads pre-assembled in the Soft Shoe. For treatment of an area other than the legs and feet, remove the Therapy Pads by sliding them through the slit on the side of the Soft Shoe. **To re-assemble the Therapy Pads into the Soft Shoe for treatment of the legs and feet follow the steps below:**

Step 1: To ensure effective treatment, you may need to adjust the placement of the Therapy Pads in the Soft Shoe to fit your foot, as shown in the images below. One Therapy Pad should be placed horizontally across the toes, and the second Therapy Pad should be placed lengthwise in the arch of your foot.

For larger feet, you may need to move the Therapy Pads away from each other.



For smaller feet, you may need to move the Therapy Pads closer to each other.



Step 2: It is very important to ensure that the diodes are making contact with your skin. Therefore, ensure that the Therapy Pad in the arch of the foot is bent toward your foot, as shown below.



INSTRUCTIONS FOR USE

Step 3: Lightly rest your foot inside the Soft Shoe. Lift the sides of the Soft Shoe and attach the top Velcro strap to the opposite side of the shoe. Then, attach the bottom Velcro strap to the opposite side of the shoe, and lastly wrap the back Velcro strap around the back of your ankle and attach the strap to the other side of the Soft Shoe.

Step 4: Never stand on the Soft Shoe. Rest your foot in the Soft Shoe or put your foot up, being sure the entire area of both Therapy Pads is touching your foot properly. Improper placement of the Therapy Pads could result in therapy failure.

Step 5: Turn on your Model 120 by depressing the ON/OFF switch and begin your treatment following the times listed in the Instructions for Use on pages 12 and 13 of this guide.

SETTING UP YOUR LEG ATTACHMENT STRAP

The Leg Attachment Strap helps facilitate the placement of the Therapy Pads around the calf or leg.

Step 1: Velcro the Therapy Pads to the Leg Attachment Strap. To ensure effective treatment, the Therapy Pads should be placed on each side of the calf, wherever the pain or poor circulation begins (as demonstrated in the picture to the right). Avoid bony areas like the shin and ankle. You may need to adjust the location of the Therapy Pads on the Leg Attachment Strap to fit your leg accordingly.

Step 2: Secure the attachment strap around your leg. Be sure not to apply the strap too tightly as this will increase your risk for a burn.



CLEANING INSTRUCTIONS

Attachment Straps:

1. Cleaning is necessary only if the Attachment Straps become excessively soiled.
2. Remove the Attachment Straps from the Therapy Pads by placing one hand on the Strap near the Velcro and gently pulling the Therapy Pad away from the strap. Note the positioning of the Attachment Straps to the Therapy Pads so you can correctly re-attach them after cleaning.
3. Hand wash the Attachment Straps in cold water and mild detergent, rinse and air dry.
4. Re-attach the Attachment Straps to the Therapy Pads by pressing the Velcro on the Attachment Straps to the Therapy Pads.

Therapy Pads

1. Cleaning is necessary only if the Therapy Pads become excessively soiled.
2. Unplug the Model 120 from the electrical outlet.
3. Lightly spray a clean cloth with a common household disinfectant spray and wipe the cloth over the entire surface of the Therapy Pads.

Soft Shoe

1. Cleaning is necessary only if the Soft Shoe becomes excessively soiled.
2. Remove the Therapy Pads from the Soft Shoe.
3. Engage the Velcro straps of the Soft Shoe.
4. Hand wash in cold or warm water with a mild detergent, rinse and air dry.
5. Take care to keep the Velcro areas free from lint.
6. Re-attach the Therapy Pads to the soft shoe as shown on page 9 of this user manual.



Do not saturate or soak the Therapy Pads in a cleaning or disinfection solution. This may cause corrosion to the Therapy Pads, which may increase the risk of burns.

PAD PLACEMENT EXAMPLES



The Therapy Pads can be arranged to fit any area of the body where the user is experiencing pain or poor circulation.



If there is an open wound in the treatment area, or if this system is being shared by more than one user;

Cover the Therapy Pads with a clear plastic barrier to avoid the spread of bacteria, viruses and other microorganisms from one person to another.

TREATMENT GUIDELINES


- Be sure you review all the **Warnings and Precautions** before using this system to prevent burns or electrical shocks.
- You may treat an area up to three times per day, separated by at least four hours between treatments.
- Treated area may be pink after therapy session, and wounds in or near the treatment area may have increased exudate.
- It is normal to experience tingling and pulsing may be felt as circulation returns to this area. If the treatment becomes uncomfortable, discontinue use and call 1-800-521-6664 for assistance.
- Diabetic patients, who may be increasing activity, should monitor their blood sugars regularly.
- Ensure that there is nothing between the Therapy Pads and your skin except for a clear plastic barrier.
 - Ensure there are no socks or any other type of clothing between the Therapy Pads and the patient's skin.
 - Do not place the Therapy Pads over a bandage, dressing or cast.
- Ensure all of the diodes on the Therapy Pad are in light, direct contact with the skin.
- This system may be safely used over metal implants, pins and screws and in people with pacemakers and defibrillators.
- It is normal to feel warmth during a therapy session. If any of the Therapy Pads feel too hot or uncomfortable, discontinue use and call 1-800-521-6664 for assistance.

Always follow your doctor's recommendations.



TREATMENT GUIDELINES

Treatment of impaired circulation, pain, stiffness or muscle spasm in the extremities

Therapy Pad Placement	Feet and Legs: - Place two Therapy Pads in the form of a "T" inside the Soft Shoe. - Place one pad on each side of calf approximately where symptoms begin using the Leg Attachment Strap.	 <p><i>Pads should be placed inside the Soft Shoe like a "T" and placed on each side of the calf.</i></p>
	Arms and Hands: - Place one to two Therapy Pads on the palm of the hand. - Place one pad on each side of the forearm or approximately where symptoms begin.	
	<p><i>Do not apply the Therapy Pads with excessive pressure</i></p>	
	Initial Usage Frequency (First 90 days)	Daily for three months or until condition resolves per a physician. <i>Using the system less than three times per week may result in therapy failure.</i>
Maintenance Frequency (Post 90 days)	1-3 times per week after initial 90-day period.	
Treatment Time	30-40 minutes	

TREATMENT GUIDELINES

Treatment of acute or chronic painful soft tissue injuries, muscle spasm or stiffness

Therapy Pad Placement	Place a minimum of two Therapy Pads directly over the painful area (i.e., knee, elbow, neck, back, etc.) If placement of Therapy Pads causes increased pain, place pads right next (lateral) to the injured site. <p><i>Do not apply the Therapy Pads with excessive pressure</i></p>
Initial Usage Frequency	Daily if possible. Minimum of three times per week.
Maintenance Usage Frequency	<i>Acute:</i> treatment may be ceased when condition resolves. <i>Chronic:</i> as needed for pain relief and circulatory maintenance.
Treatment Time	20 - 45 minutes <i>30 minutes for patients with poor skin integrity or an open wound in/near the treatment area.</i>

Therapy Pads may be easily removed from the Soft Shoe for treatment of another area of the body where needed. For Therapy Pad placement examples, reference page 14 of this user manual.